# **Low and Mid Flow Nitrogen Generators**


- Recommended and used by all major LC/MS manufacturers
- ▲ Eliminates the need for costly, dangerous, inconvenient nitrogen cylinders in the laboratory
- ▲ Models N2-04, N2-14, N2-22, N2-35 require no electricity
- Compact design frees up valuable laboratory floor space
- Phthalate-free, no organic vapors
- Unlike PSA technology, membrane will not suppress corona needle discharge.



Model N2-14 Low Flow Membrane Nitrogen Generator

Parker Balston® Low Flow Nitrogen Generators include models N2-04, N2-14, N2-14A that produce up to 61 SLPM of compressed nitrogen, on-site. The Parker Balston® Mid-Flow Nitrogen Generators include models N2-22, N2-22ANA, N2-35, and N2-35ANA that produce 132 SLPM of compressed nitrogen, on-site. The purity level of the nitrogen stream is defined by the user, for the application, and may range from 95% to 99.5%.

Low Flow Model N2-14ANA and Mid Flow Models N2-22ANA and N2-35ANA Nitrogen Generators include an oxygen analyzer which monitors the oxygen concentration of the nitrogen stream. An audible alarm signals high or low oxygen concentrations. Parker Balston Nitrogen Generators are complete systems engineered to transform standard com-



pressed air into nitrogen at safe, regulated pressures, on demand, without the need for operator attention. The systems eliminate the need for costly, dangerous dewars and cylinders in the laboratory.

Nitrogen is produced by utilizing a combination of filtration and membrane separation technologies. A high efficiency

prefiltration system pretreats the compressed air to remove all contaminants down to 0.01 micron. Hollow fiber membranes subsequently separate the clean air into a concentrated nitrogen output stream and an oxygen enriched permeate stream, which is vented from the system. The combination of these technologies produces a continuous on demand supply of pure nitrogen.



Typical applications include: LC/MS, nebulizer gas, chemical and solvent evaporation, instrument purge and supply, evaporative light scattering detector use (ELSD), and

sparging.

# **Low and Mid Flow Nitrogen Generators**

### **Nitrogen Purity / Flow Chart**

Flow measured in SLPM at indicated Operating Pressure, psig. Flows for Model N2-04 printed in black, flows for Models N2-14 and N2-14A in red.

	1	45	1	25	1	10	1	00	9	0	;	80	70	)	6	0
99.5	_	11	-	10	-	9	-	8	-	7	-	6	-	5	-	4
99	6	18	5	16	5	15	4	13	4	11	3	10	3	8	2	7
98	11	29	10	25	9	25	8	20	7	18	6	16	5	13	4	11
97	15	40	13	34	13	33	10	27	9	25	8	21	7	18	6	15
96	20	50	17	43	16	42	13	34	12	31	10	26	9	22	7	19
95	24	60	21	52	20	51	17	42	15	37	13	32	11	28	9	24

### **Nitrogen Purity / Flow Chart**

Flow measured in SLPM at indicated Operating Pressure, psig. Flows for Model N2-22, N2-22A printed in black, flows for Models N2-35, N2-35A in red.

	1	145	1	25	1	10	1	00	9	0	8	80	70		60	0
99.5	19	29	16	25	14	22	13	20	12	18	10	16	9	13	17	11
99	29	44	25	37	22	33	20	30	18	27	15	23	13	20	11	17
98	44	66	38	57	34	51	30	46	27	41	24	36	20	30	17	26
97	59	83	50	74	45	65	40	57	36	52	31	46	26	40	23	35
96	73	109	63	94	56	84	50	75	45	67	39	59	32	50	27	42
95	88	131	177	114	69	102	61	90	55	81	48	71	41	60	35	52

#### **Principal Specifications**

Model

N2-04, N2-14, N2-14ANA, N2-22, N2-22ANA, N2-35 and N2-35ANA

Nitrogen Purity 95.0% - 99.5%

 $\begin{array}{ll} \mbox{Atmospheric Dewpoint} & -58\mbox{°F (-50\mbox{°C})} \\ \mbox{Suspended Liquids} & \mbox{None} \\ \mbox{Particles} > 0.01 \mu m & \mbox{None} \\ \end{array}$ 

Commercially Sterile Yes
Phthalate-free Yes
Hydrocarbon-free Yes

Min./Max. Operating Pressure 60/145 psig

 $\begin{array}{ll} \text{Max. Press. Drop @ 99\% N}_2 \text{ Purity, 125 psig} & 10 \text{ psig} \\ \text{Recommended Ambient Operating Temperature} & 68°F (20°C) \\ \text{Max. Inlet Air Temperature} & 110°F (43°C) \\ \end{array}$ 

Inlet/Outlet Ports 1/4" NPT

Electrical Requirements N2-04, N2-14, N2-22, N2-35 None

N2-14ANA, N2-22ANA, N2-35ANA 120 VAC/60 Hz/25 Watts

Shipping Weight N2-04 42.5 lbs (19 kg)

N2-14 75 lbs (34 kg) N2-14ANA, N2-22, N2-22ANA 80 lbs (36 kg)

N2-35, N2-35ANA 90 lbs (41 kg)
Included with Model N2-14ANA, N2-22ANA, N2-35ANA

 Oxygen Analyzer
 Included with Model N2-14ANA, N2-22ANA, N2-35ANA

 Dimensions, N2-04
 16.1"h x 10.7"w x 13.4"d (40.9cm x 27.2cm x 34cm)

 Dimensions, N2-14, N2-14ANA, N2-22, N2-22ANA, N2-35, N2-35ANA
 51.5"h x 18"w x 16.2"d (130.8cm x 45.7cm x 41.1cm)

#### Ordering Information for assistance, call 800-343-4048, 8 to 5 Eastern Standard Time

Description	Galvanic Cell	Annual Maintenance Kit	Installation Kit	Preventative Maintenance Contract	Extended Support with 24 Month Warranty
N2-04	N/A	MK7840	IK7572	MFZATOC -PM	N2-04-DN2
N2-14	N/A	MK7572C	IK7572	LFMEMN2-PM	N2-14-DN2
N2-14ANA	72695A	MK7572C	IK7572	LFMEMN202-PM	N2-14ANA-DN2
N2-22, N2-35	N/A	MK7572C	IK7572	LFMEMN2-PM	N2-22-DN2, N2-35-DN2
N2-22ANA, N2-35ANA	72695A	MK7572C	IK7572	LFMEMN202-PM	N2-22ANA-DN2, N2-35ANA-DN2



## **High Flow Nitrogen Generators**


- Recommended and used by all major LC/MS manufacturers
- ▲ Eliminates the need for costly, dangerous, inconvenient nitrogen cylinders in the laboratory
- ▲ Models N2-45, N2-80, and N2-135 require no electricity
- ▲ Compact design frees up valuable laboratory floor space
- Phthalate-free, no organic vapors
- Unlike PSA technology, membrane will not suppress corona needle discharge.



Model N2-135 High Flow Membrane Nitrogen Generator

**Parker Balston® High Flow Nitrogen Generators** include models N2-45, N2-80, N2-135 that produce up to 467 SLPM of compressed nitrogen, on-site. The purity level of the nitrogen stream is defined by the user, for the application, and may range from 95% to 99.5%.

High Flow Model N2-45ANA, N2-80ANA, and N2-135ANA Nitrogen Generators include an oxygen analyzer which monitors the oxygen concentration of the nitrogen stream. An audible alarm signals high or low oxygen concentrations. Parker Balston Nitrogen Generators are complete systems engineered to transform standard compressed air into nitrogen at safe, regulated pressures, on demand, without the need for opera-



tor attention. The systems eliminate the need for costly, dangerous dewars and cylinders in the laboratory.

Nitrogen is produced by utilizing a combination of filtration and membrane separation technologies. A high efficiency prefiltration system pretreats the compressed air to remove all contaminants down to 0.01 micron. Hollow fiber membranes subsequently separate the clean air into a concentrated nitrogen output stream and an oxygen enriched permeate stream, which is vented from the system. The combination of these technologies produces a continuous on demand supply of pure nitrogen.

Typical applications include: LC/MS, nebulizer gas, chemical and solvent evaporation, instrument purge and supply, evaporative light scattering detector use (ELSD), and sparging.

# **High Flow Nitrogen Generators**

#### **Nitrogen Purity / Flow Chart** Flow LPM (liters per minute), at 68°F (25°C) inlet air temperature and operating pressure, PSIG. Flows printed in black are for Models N2-45 and N2-45A Flows printed in red are for Models N2-80 and N2-80A Flows printed in green are for Models N2-135 and N2-135A 59 78 99.5 67 100 133 83 110 92 138 183 79 106 129 194 258 73 110 147 93 124 75 101 163 244 325 94 141 187 79 119 159 97 130 200 300 400 114 171 228 97 145 194 80 119 159 233 350 467 134 201 268 111 167 222 90 135 180

Model		N2-45, N2-80, N2-135, N2-45ANA, N2-80ANA, and N2-135ANA					
Nitrogen Purity		95.0% - 99.5%					
Atmospheric Dewpoint		-58°F (-50°C)					
Suspended Liquids		None					
Particles > 0.01µm		None					
Commercially Sterile		Yes					
Phthalate-free		Yes					
Hydrocarbon-free		Yes					
Min./Max. Operating Pre	essure	60/145 psig					
Max. Press. Drop @ 99	% N <sub>2</sub> Purity, 125 psig	10 psig					
Recommended Ambient	Operating Temperature	72°F (22°C)					
Max. Inlet Air Temperatu	ire	110°F (43°C)					
Inlet/Outlet Ports		1/2" NPT					
<b>Electrical Requirements</b>	N2-45, N2-80, N2-135	None					
	N2-45ANA, N2-80ANA, N2-135ANA	120 VAC/60 Hz/25 Watts					
Shipping Weight	N2-45, N2-80, N2-135	250 lbs (114 kg)					
	N2-45ANA, N2-80ANA, N2-135ANA	250 lbs (114 kg)					
Oxygen Analyzer		Included with Model N2-45ANA, N2-80ANA, N2-135ANA					
Dimensions		67"h x 24"w x 20"d (140cm x 61cm x 50cm)					

Ordering Information for assistance, call 800-343-4048, 8 to 5 Eastern Time									
Description	Galvanic Cell	Carbon Tower	Installation Kit	Preventative Maintenance Contract	Extended Support with 24 Month Warranty				
N2-45	N/A	75344	IK75880	MFMEMN2-PM	N2-45-DN2				
N2-45ANA	72695A	75344	IK75880	MFMEMN202-PM	N2-45ANA-DN2				
N2-80	N/A	75344	IK75880	MFMEMN2-PM	N2-80-DN2				
N2-80ANA	72695A	75344	IK75880	MFMEMN202-PM	N2-80ANA-DN2				
N2-135	N/A	75344	IK75880	MFMEMN2-PM	N2-135-DN2				
N2-135ANA	72695A	75344	IK75880	MFMEMN202-PM	N2-135ANA-DN2				

